

# The C-MAC<sup>®</sup> Pocket Monitor

Portability with real-time documentation



# C-MAC® Pocket Monitor

## The premium class of video laryngoscopy

The C-MAC® POCKET MONITOR sets new standards in image quality for mobile video laryngoscopy. Its unique display allows intuitive work, even in emergency situations.

- 3.5" OTI® display with wide viewing angle of 160°
- Real-time documentation of images and videos via the BlueButton
- Equipped with the C-MAC® system interface, ready for future applications
- IPX8: Absolutely water- and dustproof, validated for wipe disinfection



### C-MAC® video laryngoscopes (reusable)

- Standard MACINTOSH and MILLER blades for safe intubation of infants, children, and adults
- Highly curved D-BLADE and D-BLADE Ped for the particularly difficult airway
- The highest ergonomic requirements – specially designed and shortened handle
- Robust and lightweight: High stability and resistance to impact; light metal design



### C-MAC® S video laryngoscopes (for single use)

- C-MAC® S IMAGER with C-MAC® system interface
- C-MAC® S video laryngoscopes available in adult sizes in MAC #3, MAC #4 and D-BLADE versions and Pediatric and Neonate sizes in MILLER 0 and 1 versions
- Optimal protection against cross-contamination as blade and handle form one continuous piece
- Extremely high tensile strength
- Compact design with ergonomic handle



### Premium accessory

- The C-MAC® PM connecting cable enables flexible positioning of the monitor for difficult-to-access patients
- Uninterrupted readiness at no extra cost thanks to exchangeable and rechargeable battery
- Foreign body forceps specially adapted to the MACINTOSH blade and curved video laryngoscopes
- Special bag designed for the C-MAC® POCKET MONITOR, C-MAC® video laryngoscopes and accessories



# C-MAC<sup>®</sup> PM – can be used at any time and in any place



Intensive Care Medicine



Prehospital Care



General Anesthesia



Emergency Room

It is recommended to check the suitability of the product for the intended procedure prior to use.





...evolution continues

**STORZ**  
**KARL STORZ — ENDOSKOPE**

*THE DIAMOND STANDARD*

KARL STORZ SE & Co. KG  
Dr.-Karl-Storz-Straße 34, 78532 Tuttlingen/Germany  
Postbox 230, 78503 Tuttlingen/Germany  
Phone: +49 (0)7461 708-0  
Fax: +49 (0)7461 708-105  
E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
[www.karlstorz.com](http://www.karlstorz.com)